

OTOair - wireless made easy



Increased ease of use and convenience for your wireless Otometrics instruments

- Plug and play operation
- Independent of other Bluetooth installations
- Supported by OTOSuite version 3.50 or higher

OTOair is a USB dongle for OTOSuite compatible Bluetooth instruments, i.e. MADSEN OTOflex 100, AURICAL SpeechLink 100, AURICAL VisibleSpeech and AURICAL FreeFit.

Plug and Play

The OTOair dongle plugs directly into any compatible USB port on your computer for fast and easy operation making it a great solution for working between different work-stations. You don't even have to remember to insert it into the same USB port – your computer will recognize it every time. Scanning for devices is faster with OTOair as only Otometrics devices are recognized in the scanning process.

“At one customer installation it was not possible to get a connection to any Bluetooth device. But with the OTOair we were able to make a pairing to the FreeFit and the customer could finally start working” states an Otometrics service technician in Germany.

OTOair does not interfere with or depend on Windows Bluetooth stack or 3rd party drivers, which eliminates the problems and troubleshooting known from traditional Bluetooth setups.

New orders of the listed devices will automatically include an OTOair – existing customers that would like to experience the benefits of OTOair can purchase it as an accessory (part no. 8-62-41503).