



Examination Table

Model ET-600



Specifically Designed for VNG and ENG Testing

- Comfortable
- Convenient
- Economical
- Attractive



Specifically Designed for ENG Testing

Examination Table

Model ET-600

KEY FEATURES

- Headrest for comfortable patient positioning, removable for Hallpike maneuvers or head hanging positions.
- Precise 30° positioning for accurate caloric testing and easy access to ears during irrigations.
- Comfortable 2-inch foam padding covered by leatherette.
- Attractive dark wood base.
- Large storage capacity beneath the table. (Fully enclosed base is also available, Model ET-600A or fully enclosed base with shelf, Model ET-600AS.)

SPECIFICATIONS

Size (W x L x H)

30 x 78 x 30 inches

76 x 198 x 76 centimeters

Weight

148lbs.

67kg.



INNOVATION IN PRACTICE

ICS Medical Corporation
125 Commerce Drive
Schaumburg, IL 60173-5329

tel 800 289 2150
tel 847 534 2150
fax 847 534 2151
sales@icsmedical.com
www.icsmedical.com