

GN Otometrics

Product news

13.05.04

OTOdiagnostics Suite - new audiologic software from GN Otometrics

GN Otometrics recently introduced OTOdiagnostics Suite, versatile and intuitive audiologic software for interfacing with our new generation of hearing test devices. Currently, the Suite is supplied together with the OTOflex 100, a compact wireless middle-ear analyser capable of performing all diagnostic and screening immittance tests.

The software represents an exciting new conceptual framework for test and measurement software. It permits remote operation of the OTOflex 100 from a PC and printout of all test results on a single sheet of paper. If the middle-ear analyser is used in the proximity of a PC, test results will be displayed on the PC's screen at the same time as on the analyser's display. If the analyser is used elsewhere in the clinic or hospital, measurement data can be transferred to the PC or printer *without the two devices having to be physically connected*. All data transfer is wireless, via Bluetooth.

OTOdiagnostics Suite can be run stand-alone or as a measurement module from a NOAH 3 compatible practice management system such as NOAH System 3.1 or GN Otometrics' Mirage or AuditBase Systems. Access to a NOAH 3 compatible database also turns OTOdiagnostics Suite into a viewer for immittance test results acquired from the MADSEN Zodiac.

Readers may contact their local GN Otometrics' representative for a demonstration of the OTOdiagnostics Suite.

Contact: GN Otometrics

Tel. +45 72 111 555

Email: info@gnotometrics.dk

Website: www.gnotometrics.com

