

FOR IMMEDIATE RELEASE

Contact: GN Otometrics Sales, North America

sales@gnotometrics.com

www.gnotometrics.com

1 800-289-2150

SCHAUMBURG, Ill., 10 June 2005

ICS CHARTR EP w/ASSR sw-upgrade

The upgrade package (ver. 5.2) includes new enhancements to our Auditory Steady State Response. In our continuing effort to express our competitiveness in this market new features has been added. This includes features such as 2 new frequencies (250Hz and 8 kHz), user modification of ASSR Correction Factors, user defined ASSR protocols and Masking. Also, EP features include incorporated ability to review databases archived on CD-ROM and enhanced report layout.